

CONTENTS LISTING

TMJ Fixation Instrument Set (Loaner) Part Number 60-0020

Please return set in a clean and sterilized condition to:
TMJ Concepts
 6059 Knoll Drive
 Ventura, CA 93003 USA

This set is the property of TMJ Concepts. No return authorization is needed. Retain all unused bone screws in set. Discard any screws that have had patient contact during surgery.

INTERNATIONAL INSTRUMENT SET RETURN INFORMATION:

TMJ Concepts is the original manufacturer of this instrument set. The contents and materials are described below. The complete set is valued at \$4500 USD, however it has been loaned free of charge to a doctor or hospital for use in implanting TMJ implant components.

US FDA Facility Registration Number: 2031049
 US FDA Medical Device Listing Number: MDL #A962379
 US FDA Premarket Approval Number: PMA P980052

INSTRUMENTS

Qty	Part No.	Description	Material
2	60-0100	Quick-Connect Instrument Handle	SS / Al / Acetal Copolymer
2	60-0110	Driver Blade	Stainless Steel
1	60-0120	Implant Stabilizer	Stainless Steel
1	60-0130	Mandibular Forceps	Stainless Steel
1	60-0200	Fossa Seating Tool	Stainless Steel / Co-Cr
1	60-0300	Drill Guide	Stainless Steel
2	60-0420	Pilot Drill 2.0mm	Stainless Steel
2	60-0423	Pilot Drill 2.3mm	Stainless Steel
1	60-0500	Sterilization Case Lid	Polyphenylsulfone (PPSU)
1	60-0510	Sterilization Case Upper Tray	Polypropylene / PPSU / SS
1	60-0520	Sterilization Case Lower Tray	Polypropylene / SS

BONE SCREWS

Qty*	Part No.	2.0mm Diameter Screws	Material
0	61-2005	Bone Screw 2.0mm x 5mm	Titanium 6Al-4V Alloy
21	61-2006	Bone Screw 2.0mm x 6mm	Titanium 6Al-4V Alloy
21	61-2008	Bone Screw 2.0mm x 8mm	Titanium 6Al-4V Alloy
28	61-2010	Bone Screw 2.0mm x 10mm	Titanium 6Al-4V Alloy
14	61-2012	Bone Screw 2.0mm x 12mm	Titanium 6Al-4V Alloy
7	61-2014	Bone Screw 2.0mm x 14mm	Titanium 6Al-4V Alloy
7	61-2016	Bone Screw 2.0mm x 16mm	Titanium 6Al-4V Alloy
Qty*	Part No.	2.3mm Diameter Safety Screws	Material
6	61-2306	Bone Screw 2.3mm x 6mm	Titanium 6Al-4V Alloy
6	61-2308	Bone Screw 2.3mm x 8mm	Titanium 6Al-4V Alloy
12	61-2310	Bone Screw 2.3mm x 10mm	Titanium 6Al-4V Alloy
6	61-2312	Bone Screw 2.3mm x 12mm	Titanium 6Al-4V Alloy
6	61-2314	Bone Screw 2.3mm x 14mm	Titanium 6Al-4V Alloy
6	61-2316	Bone Screw 2.3mm x 16mm	Titanium 6Al-4V Alloy

* Minimum initial screw quantity (Additional screws may be present in set.)

Bone screws are made from implant grade Ti-6Al-4V ELI (ASTM F136)

INSTRUMENT SET LAYOUT

